



NUMBER 1 IN COMPRESSION THERAPY

SPECIALTY MEDICAL PRODUCTS, INC.

Billing: 50 Pennsylvania Ave Malvern, PA 19382
Shipping: 187 Pennsylvania Ave Malvern, PA 19382
Phone: (888) 698-1522 Fax (610) 644-3992
Email: info@smplecares.com
Website: www.smplecares.com

Sequential Circulator Model SC-2004



INTENDED USE Model SC-2004 Sequential Circulator is a gradient, sequential, pneumatic compression device, intended for the primary or adjunctive treatment of primary or secondary Lymphedema. The device is intended for alternate treatment of chronic venous stasis ulcers and associated venous insufficiency, as well as general treatment of swelling of the extremities. The device is intended for home or hospital use.

Specifications:

- Dimensions: 5.5 x 12 x 8 (H x W x D)
- Weight: 8 lbs
- Pressure Range: 0-125mm Hg
- Electrical Rating: 120 VAC, 60 HZ, .5 A
- Cycle Time: 18 Seconds/ Chamber
- Inflation: 72 secs
- Deflation: 18 secs
- Warranty: 3 years
- Bi-lateral operation:



Before



After



What Does it Do?

1. Reduces Pain
2. Reduces Edema
3. Mimics the Lymphatic System
4. Promotes Lymphatic Flow by Moving Fluid in the Proper Physiological Direction (distal to proximal)
5. Stimulates Circulation to Promote Wound Healing
6. Delivers Continuous Flow
7. Bilateral Operation
8. Lightweight . Portable . Easy to Use . Cost-Effective

CONTRAINDICATIONS

Pneumatic compression is contraindicated for the patients with:

1. Congestive heart failure
2. Deep Vein Thrombosis
3. Inflammatory phlebitis or episodes of pulmonary embolism
4. Infections in the limb, including cellulitis, without appropriate antibiotic coverage
5. Presence of lymphangiosarcoma

Bio Compression Systems

Lower Extremity/ 4-Chamber Garment Sizing Chart

Part Number	Size	Circumferences							Lengths		
		Thigh below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease		
GS-3045-H	Half	Standard	NA	29"	27"	23.5"	22"	19"	NA		
GW-3045-H		Wide		30.5"	28"	29"	23"				
GA-3045-H		Adjustable		32"	30"	23.5"	22"				
GWA-3045-H		Adj. Wide		33.5"	31"	29"	23"				
GN-3045-S	Small	Narrow	25.5"	25"	22"	20"	19.5"	25.5"	27"		
GS-3045-S		Standard	31"	30"	28"	24.5"	19.5"				
GW-3045-S		Wide	35.5"	34"	30.5"	28"	21.5"			21.5"	
GA-3045-S		Adjustable	34"	33"	31"	27.5"	19.5"			19.5"	
GWA-3045-S		Adj. Wide	38.5"	37"	33.5"	31"	25.5"	21.5"			
GN-3045-M	Medium	Narrow	27.5"	25.5"	23"	20"	19.5"	28"	32"		
GS-3045-M		Standard	32"	30"	28"	24.5"	19.5"				
GW-3045-M		Wide	36.5"	34"	30"	28"	22.5"			22.5"	
GXA-3045-M		Extra Wide	44"	41"	36"	34"	30"			22.5"	22.5"
GA-3045-M		Adjustable	35"	34"	31"	27.5"	19.5"	19.5"			
GWA-3045-M		Adj. Wide	39.5"	36.5"	33.5"	31"	22.5"	22.5"			
GXA-3045-M		Adj. Extra Wide	47"	44"	39"	37"	22.5"	22.5"			
GN-3045-L	Large	Narrow	28"	26"	24"	21"	19.5"	32"	36"		
GS-3045-L		Standard	33"	32"	30"	26"	19.5"				
GW-3045-L		Wide	37"	36"	32"	29"	20.5"			20.5"	
GA-3045-L		Adjustable	37"	36"	33"	29"	19.5"			19.5"	
GWA-3045-L		Adj. Wide	40"	38"	35"	32"	22.5"	22.5"			